



TTC23 Monoclonal Antibody

Catalog No	YP-mAb-04261
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TTC23
Immunogen	The antiserum was produced against synthesized peptide derived from human TTC23. AA range:341-390
Specificity	TTC23 Monoclonal Antibody detects endogenous levels of TTC23 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cervical cancer proto oncogene 8 protein;Cervical cancer proto-oncogene 8 protein;HCC 8;HCC-8;OTTHUMP00000190851;Protooncogene 8;TTC23;Tetratricopeptide repeat protein 23;HCC8;TPR repeat protein 23
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	Cell projection, cilium . Colocalizes with EVC and IQCE at the EvC zone of primary cilia. .
Tissue Specificity	Colon,Coronary artery,Placenta,
Function	similarity:Contains 4 TPR repeats.,
Background	similarity:Contains 4 TPR repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

