



β -1,4-Gal-T5 Monoclonal Antibody

Catalog No	YP-mAb-04294
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	B4GALT5
Immunogen	The antiserum was produced against synthesized peptide derived from human B4GALT5. AA range:321-370
Specificity	β-1,4-Gal-T5 Monoclonal Antibody detects endogenous levels of β-1,4-Gal-T5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	4-galactosyltransferase 5;4-GalT II;4-GalTase 5;B4Gal T5;b4Gal-T5;B4galt5; B4GT5_HUMAN;Beta 1 4 galactosyltransferase 5;Beta 1 4 galactosyltransferase V;Beta 1 4 GalT II;Beta 1 4 GalT IV;Beta 1 4 GalTase 5; Beta-1;Beta4Gal-T5;UDP-Gal:beta-GlcNAc beta-1;UDP-galactose:beta-N-acetylglucosamine beta-1
Observed Band	40kD
Calculated Molecular Weight	40kD
Cell Pathway	Golgi apparatus, Golgi stack membrane ; Single-pass type II membrane protein. Golgi apparatus . Trans cisternae of Golgi stack. .
Tissue Specificity	Ubiquitously expressed.
Function	cofactor:Manganese.,function:Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids.,online information:Beta-1,4-galactosyltransferase 5, online information:GlycoGene database,pathway:Protein modification; protein glycosylation.,similarity:Belongs to the glycosyltransferase 7 family.,subcellular location:Trans cisternae of Golgi stack.,tissue specificity:Ubiquitously expressed.,



Background

This gene is one of seven beta-1,4-galactosyltransferase (beta4GalT) genes. They encode type II membrane-bound glycoproteins that appear to have exclusive specificity for the donor substrate UDP-galactose; all transfer galactose in a beta1,4 linkage to similar acceptor sugars: GlcNAc, Glc, and Xyl. Each beta4GalT has a distinct function in the biosynthesis of different glycoconjugates and saccharide structures. As type II membrane proteins, they have an N-terminal hydrophobic signal sequence that directs the protein to the Golgi apparatus and which then remains uncleaved to function as a transmembrane anchor. By sequence similarity, the beta4GalTs form four groups: beta4GalT1 and beta4GalT2, beta4GalT3 and beta4GalT4, beta4GalT5 and beta4GalT6, and beta4GalT7. The function of the enzyme encoded by this gene is not clear. This gene was previously designated as B4GALT4 but was renamed


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using β -1,4-Gal-T5 Monoclonal Antibody