



KIF1C Monoclonal Antibody

Catalog No	YP-mAb-04889
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIF1C KIAA0706
Immunogen	Synthesized peptide derived from human protein . at AA range: 1110-1190
Specificity	KIF1C Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0706;KIF1C;kinesin family member 1C;Kinesin like protein KIF1C;LTXS1;KIF1C KIAA0706;Kinesin-like protein KIF1C
Observed Band	121kD
Calculated Molecular Weight	121kD
Cell Pathway	Cytoplasm, cytoskeleton .
Tissue Specificity	Expressed in all tissues examined, with most abundant expression in heart and skeletal muscle.
Function	function:Motor required for the retrograde transport of Golgi vesicles to the endoplasmic reticulum. Has a microtubule plus end-directed motility., PTM:Phosphorylated on tyrosine residues.,similarity:Belongs to the kinesin-like protein family. Unc-104 subfamily.,similarity:Contains 1 FHA domain.,similarity:Contains 1 kinesin-motor domain.,subunit:Monomer .,tissue specificity:Expressed in all tissues examined, with most abundant expression in heart and skeletal muscle.,
Background	kinesin family member 1C(KIF1C) Homo sapiens The protein encoded by this gene is a member of the kinesin-like protein family. The family members are microtubule-dependent molecular motors that transport organelles within cells



and move chromosomes during cell division. Mutations in this gene are a cause of spastic ataxia 2, autosomal recessive. [provided by RefSeq, May 2014],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using KIF1C Monoclonal Antibody