



KKCC1 Monoclonal Antibody

Catalog No	YP-mAb-04890
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CAMKK1 CAMKKA
Immunogen	Synthesized peptide derived from human protein . at AA range: 400-480
Specificity	KKCC1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CaM KK 1;CaM kinase kinase alpha;CaM kinase kinase 1;CaM kinase IV kinase; Calcium/calmodulin-;KKCC1;CAMKK1 CAMKKA;Calcium/calmodulin-dependent protein kinase kinase 1 (CaM-KK 1) (CaM-kinase kinase 1) (CaMKK 1) (EC 2.7.11.17) (CaM-kinase IV kinase) (Calcium/calmodulin-dependent protein kinase kinase alpha) (CaM-KK
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Amygdala,Brain,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,domain:The autoinhibitory domain overlaps with the calmodulin binding region and may be involved in intrasteric autoinhibition.,domain:The RP domain (arginine/proline-rich) is involved in the recognition of CAMKI and CAMK4 as substrates.,enzyme regulation:Activated by Ca(2+)/calmodulin. Binding of calmodulin may release intrasteric autoinhibition. Partially inhibited upon phosphorylation by PRCAPA/ PKA (By similarity). May be regulated through phosphorylation by CAMK1 and CAMK4.,function:Calcium/calmodulin-dependent protein kinase that belongs to a proposed calcium-triggered signaling cascade involved in a number of cellular



processes. Phosphorylates CAMK1, CAMK1D, CAMK1G and CAMK4. Involved in regulating cell apoptosis. Promotes cell survival by phosphorylating AKT1/ PKB that inhibits pro-apoptotic BAD/Bcl2-antagonist of cell de

Background

The product of this gene belongs to the Serine/Threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. This protein plays a role in the calcium/calmodulin-dependent (CaM) kinase cascade. Three transcript variants encoding two distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using KKCC1 Monoclonal Antibody