



TF2LX Monoclonal Antibody


Catalog No	YP-mAb-04950
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TGIF2LX TGIFLX
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TF2LX Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Homeobox protein TGIF2LX;TGIF like on the X;TGIF2LX;TGIFLX;TF2LX; TGIF2LX TGIFLX;Homeobox protein TGIF2LX (TGF-beta-induced transcription factor 2-like protein) (TGFB-induced factor 2-like protein;X-linked) (TGIF-like on the X); X-linked) (TGIF-like on the X)
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Nucleus .
Tissue Specificity	Specifically expressed in adult testis.
Function	function:May have a transcription role in testis.,similarity:Belongs to the TALE/ TGIF homeobox family.,similarity:Contains 1 homeobox DNA-binding domain., tissue specificity:Specifically expressed in adult testis.,
Background	This gene encodes a member of the TALE/TGIF homeobox family of transcription factors. Testis-specific expression suggests that this gene may play a role in spermatogenesis. A homolog of this gene lies within the male specific region of chromosome Y, in a block of sequence that is thought to be the result of a large X-to-Y transposition. [provided by RefSeq, Jul 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TF2LX Monoclonal Antibody