



EFNB2 Monoclonal Antibody

Catalog No	YP-mAb-04973
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EFNB2 EPLG5 HTKL LERK5
Immunogen	Synthesized peptide derived from human protein . at AA range: 210-290
Specificity	EFNB2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EFNB2;EPLG5;EPH-related receptor tyrosine kinase ligand 5;Ephrin-B2;Htk L; EFNB2 EPLG5 HTKL LERK5;Ephrin-B2 (EPH-related receptor tyrosine kinase ligand 5) (LERK-5) (HTK ligand) (HTK-L)
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Cell junction, adherens junction .
Tissue Specificity	Lung and kidney.
Function	function: Binds to the receptor tyrosine kinases EPHB4 and EPHA3. May play a role in constraining the orientation of longitudinally projecting axons., PTM: Inducible phosphorylation of tyrosine residues in the cytoplasmic domain., similarity: Belongs to the ephrin family., subunit: Interacts with PDZRN3 (By similarity). Binds to the receptor tyrosine kinases EPHB4 and EPHA3. Binds to Nipah virus G protein., tissue specificity: Lung and kidney.,
Background	This gene encodes a member of the ephrin (EPH) family. The ephrins and EPH-related receptors comprise the largest subfamily of receptor protein-tyrosine kinases and have been implicated in mediating developmental events,



especially in the nervous system and in erythropoiesis. Based on their structures and sequence relationships, ephrins are divided into the ephrin-A (EFNA) class, which are anchored to the membrane by a glycosylphosphatidylinositol linkage, and the ephrin-B (EFNB) class, which are transmembrane proteins. This gene encodes an EFNB class ephrin which binds to the EPHB4 and EPHA3 receptors. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using EFNB2 Monoclonal Antibody