



H1X Monoclonal Antibody


Catalog No	YP-mAb-05012
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	H1FX
Immunogen	Synthesized peptide derived from human protein . at AA range: 100-180
Specificity	H1X Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Histone H1.10;Histone H1;Histone;H1X;H1-10;H1FX;Histone H1x
Observed Band	23kD
Calculated Molecular Weight	23kD
Cell Pathway	Nucleus. Chromosome.
Tissue Specificity	Expressed ubiquitously.
Function	function:Histones H1 are necessary for the condensation of nucleosome chains into higher order structures.,similarity:Belongs to the histone H1/H5 family.,tissue specificity:Expressed ubiquitously.,
Background	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene encodes a replication-independent histone that is a member of the histone H1 family. [provided by RefSeq, Oct 2015],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using H1X Monoclonal Antibody