



KPSH2 Monoclonal Antibody


Catalog No	YP-mAb-05017
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PSKH2
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	KPSH2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Serine/threonine-protein kinase H2;PSK-H2;PSK H2;Protein serine kinase H2; EC:2.7.11.1;KPSH2;PSKH2;Serine/threonine-protein kinase H2 (EC 2.7.11.1) (Protein serine kinase H2) (PSK-H2)
Observed Band	42kD
Calculated Molecular Weight	42kD
Cell Pathway	
Tissue Specificity	
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,caution:Asn-183 is present instead of the conserved Asp which is expected to be an active site residue.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,caution:Asn-183 is present instead of the conserved Asp which is expected to be an active site residue.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using KPSH2 Monoclonal Antibody