



DNJB5 Monoclonal Antibody


Catalog No	YP-mAb-05068
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNAJB5 HSC40
Immunogen	Synthesized peptide derived from human protein . at AA range: 140-220
Specificity	DNJB5 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DNJB5;DNAJB5 HSC40;DnaJ homolog subfamily B member 5 (Heat shock protein Hsp40-2) (Heat shock protein Hsp40-3) (Heat shock protein cognate 40) (Hsc40)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	nucleus,cytosol,
Tissue Specificity	Brain,Colon,Skin,
Function	induction:Expressed under normal conditions, its expression can further be increased after various stress treatments.,similarity:Contains 1 J domain.,
Background	DNAJB5 belongs to the evolutionarily conserved DNAJ/HSP40 protein family. For background information on the DNAJ family, see MIM 608375.[supplied by OMIM, Mar 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using DNJB5 Monoclonal Antibody