



FRAT1 Monoclonal Antibody


Catalog No	YP-mAb-05076
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FRAT1
Immunogen	Synthesized peptide derived from human protein . at AA range: 130-210
Specificity	FRAT1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FRAT1;Proto-oncogene FRAT1 (Frequently rearranged in advanced T-cell lymphomas 1) (FRAT-1)
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Cytoplasm.
Tissue Specificity	Lung,Testis,
Function	function:Positively regulates the Wnt signaling pathway by stabilizing beta-catenin through the association with GSK-3. May play a role in tumor progression and collaborate with PIM1 and MYC in lymphomagenesis., PTM:Phosphorylated.,similarity:Belongs to the GSK-3-binding protein family., subunit:Binds DVL1. Binds GSK-3 and prevent GSK-3-dependent phosphorylation.,
Background	The protein encoded by this gene belongs to the GSK-3-binding protein family. The protein inhibits GSK-3-mediated phosphorylation of beta-catenin and positively regulates the Wnt signaling pathway. It may function in tumor progression and in lymphomagenesis. [provided by RefSeq, Oct 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using FRAT1 Monoclonal Antibody