



GOT1B Monoclonal Antibody


Catalog No	YP-mAb-05194
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GOLT1B GOT1A CGI-141 HDCMA39P UNQ432/PRO793
Immunogen	Synthesized peptide derived from human protein . at AA range: 70-150
Specificity	GOT1B Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GOT1B;GOLT1B GOT1A CGI-141 HDCMA39P UNQ432/PRO793;Vesicle transport protein GOT1B (Golgi transport 1 homolog B) (Putative NF-kappa-B-activating protein 470) (hGOT1a)
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed. Tends to be up-regulated in seminomas compared to normal testis.
Function	function:May be involved in fusion of ER-derived transport vesicles with the Golgi complex.,similarity:Belongs to the GOT1 family.,
Background	function:May be involved in fusion of ER-derived transport vesicles with the Golgi complex.,similarity:Belongs to the GOT1 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GOT1B Monoclonal Antibody