



BSSP4 Monoclonal Antibody

Catalog No	YP-mAb-05387
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRSS22 BSSP4 PRSS26 SP001LA UNQ302/PRO343
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	BSSP4 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BSSP4;PRSS22 BSSP4 PRSS26 SP001LA UNQ302/PRO343;Brain-specific serine protease 4 (BSSP-4) (EC 3.4.21.-) (Serine protease 22) (Serine protease 26) (Trypsin epsilon)
Observed Band	34kD
Calculated Molecular Weight	34kD
Cell Pathway	Secreted .
Tissue Specificity	Expressed abundantly in the epithelial cells of the airways, including trachea, esophagus and fetal lung. Scarce in adult lung. Expressed at low levels in placenta, pancreas, prostate and thyroid gland.
Function	function:Preferentially cleaves the synthetic substrate H-D-Leu-Thr-Arg-pNA compared to tosyl-Gly-Pro-Arg-pNA.,similarity:Belongs to the peptidase S1 family.,similarity:Contains 1 peptidase S1 domain.,tissue specificity:Expressed abundantly in the epithelial cells of the airways, including trachea, esophagus and fetal lung. Scarce in adult lung. Expressed at low levels in placenta, pancreas, prostate and thyroid gland.,

**Background**


This gene encodes a member of the trypsin family of serine proteases. The enzyme is expressed in the airways in a developmentally regulated manner. The gene is part of a cluster of serine protease genes on chromosome 16. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using BSSP4 Monoclonal Antibody