



CAMLG Monoclonal Antibody

Catalog No	YP-mAb-05393
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CAMLG CAML
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CAMLG Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	calcium modulating ligand;Calcium signal-modulating cyclophilin ligand;CAML; GET2;Guided entry of tail-anchored proteins factor CAMLG;CAMLG;CAMLG CAML;Calcium signal-modulating cyclophilin ligand (CAML)
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Ubiquitous. Highest levels in brain, testis and ovary.
Function	function:Likely involved in the mobilization of calcium as a result of the TCR/CD3 complex interaction. Binds to cyclophilin B.,subunit:The N-terminal domain binds to TNFRSF13B/TACI.,tissue specificity:Ubiquitous. Highest levels in brain, testis and ovary.,
Background	The immunosuppressant drug cyclosporin A blocks a calcium-dependent signal from the T-cell receptor (TCR) that normally leads to T-cell activation. When bound to cyclophilin B, cyclosporin A binds and inactivates the key signaling intermediate calcineurin. The protein encoded by this gene functions similarly to cyclosporin A, binding to cyclophilin B and acting downstream of the TCR and upstream of calcineurin by causing an influx of calcium. This integral membrane



protein appears to be a new participant in the calcium signal transduction pathway, implicating cyclophilin B in calcium signaling, even in the absence of cyclosporin. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CAMLG Monoclonal Antibody