




SPA9 Monoclonal Antibody

Catalog No	YP-mAb-06157
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SERPINA9 GCET1 SERPINA11 UNQ692/PRO1337
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SPA9 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPA9;SERPINA9 GCET1 SERPINA11 UNQ692/PRO1337;Serp A9 (Centerin) (Germinal center B-cell-expressed transcript 1 protein)
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	[Isoform 1]: Secreted.; [Isoform 2]: Cytoplasm .; [Isoform 3]: Cytoplasm .; [Isoform 4]: Cytoplasm .; [Isoform 5]: Cytoplasm .; [Isoform 6]: Cytoplasm .; [Isoform 7]: Membrane ; Single-pass type II membrane protein .
Tissue Specificity	Highly expressed in normal germinal center (GC) B-cells and GC B-cell-derived malignancies.
Function	similarity:Belongs to the serpin family.,tissue specificity:Highly expressed in normal germinal center (GC) B-cells and GC B-cell-derived malignancies.,
Background	similarity:Belongs to the serpin family.,tissue specificity:Highly expressed in normal germinal center (GC) B-cells and GC B-cell-derived malignancies.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SPA9 Monoclonal Antibody