



ZA2G Monoclonal Antibody

Catalog No	YP-mAb-06421
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AZGP1 ZAG ZNGP1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZA2G Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AZGP1;Zinc-alpha-2-glycoprotein;Zn alpha 2 GP;Zn alpha 2 glycoprotein;ZAG;ZA2G;AZGP1 ZAG ZNGP1;Zinc-alpha-2-glycoprotein (Zn-alpha-2-GP) (Zn-alpha-2-glycoprotein)
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	Secreted.
Tissue Specificity	Blood plasma, seminal plasma, urine, saliva, sweat, epithelial cells of various human glands, liver.
Function	function:Stimulates lipid degradation in adipocytes and causes the extensive fat losses associated with some advanced cancers. May bind polyunsaturated fatty acids.,similarity:Belongs to the MHC class I family.,similarity:Contains 1 Ig-like C1-type (immunoglobulin-like) domain.,subunit:Interacts with PIP.,tissue specificity:Blood plasma, seminal plasma, urine, saliva, sweat, epithelial cells of various human glands, liver.,
Background	function:Stimulates lipid degradation in adipocytes and causes the extensive fat losses associated with some advanced cancers. May bind polyunsaturated fatty acids.,similarity:Belongs to the MHC class I family.,similarity:Contains 1 Ig-like



C1-type (immunoglobulin-like) domain.,subunit:Interacts with PIP.,tissue specificity:Blood plasma, seminal plasma, urine, saliva, sweat, epithelial cells of various human glands, liver.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ZA2G Monoclonal Antibody