



CCNY Monoclonal Antibody


Catalog No	YP-mAb-06436
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCNY C10orf9 CBCP1 CFP1
Immunogen	Synthesized peptide derived from human protein . at AA range: 270-350
Specificity	CCNY Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C10orf9;CBCP1;CFP1;Cyc Y;Cyclin box protein 1;CCNY;CCNY C10orf9 CBCP1 CFP1;Cyclin-Y (Cyc-Y) (Cyclin box protein 1) (Cyclin fold protein 1) (cyclin-X)
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side .; [Isoform 3]: Nucleus.
Tissue Specificity	Widely expressed.
Function	similarity:Belongs to the cyclin family. Cyclin Y subfamily.,similarity:Contains 1 cyclin N-terminal domain.,
Background	Cyclins, such as CCNY, control cell division cycles and regulate cyclin-dependent kinases (e.g., CDC2; MIM 116940) (Li et al., 2009 [PubMed 18060517]).[supplied by OMIM, May 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CCNY Monoclonal Antibody