




TRI31 Monoclonal Antibody

Catalog No	YP-mAb-06567
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRIM31 C6orf13 RNF
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TRI31 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TRIM;RNF;HCGI;HCG1;EC:2.3.2.27;TRI31;TRIM31 C6orf13 RNF;E3 ubiquitin-protein ligase TRIM31 (EC 6.3.2.-) (Tripartite motif-containing protein 31)
Observed Band	46kD
Calculated Molecular Weight	46kD
Cell Pathway	Cytoplasm . Mitochondrion . Predominantly expressed in the cytoplasm but a fraction is associated with the mitochondria.
Tissue Specificity	Up-regulated in gastric adenocarcinomas.
Function	similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 B box-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,
Background	This gene encodes a protein that functions as an E3 ubiquitin-protein ligase. This gene shows altered expression in certain tumors and may be a negative regulator of cell growth. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TRI31 Monoclonal Antibody