




MANs1 Monoclonal Antibody

Catalog No	YP-mAb-06575
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MANSC1 LOH12CR3 UNQ316/PRO361
Immunogen	Synthesized peptide derived from human protein . at AA range: 310-390
Specificity	MANs1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	9130403P13Rik;LOH12CR3;MANSC domain containing 1;MANSC1;MANs1;MANSC1 LOH12CR3 UNQ316/PRO361;MANSC domain-containing protein 1 (Loss of heterozygosity 12 chromosomal region 3 protein)
Observed Band	47kD
Calculated Molecular Weight	47kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Widely expressed.
Function	similarity:Contains 1 MANSC domain.,tissue specificity:Widely expressed.,
Background	similarity:Contains 1 MANSC domain.,tissue specificity:Widely expressed.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MANS1
Monoclonal Antibody