




PAR16 Monoclonal Antibody

Catalog No	YP-mAb-06803
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PARP16 C15orf30
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PAR16 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PAR16;PARP16 C15orf30;Mono [ADP-ribose] polymerase PARP16 (EC 2.4.2.31) (ADP-ribosyltransferase diphtheria toxin-like 15) (ARTD15) (Poly [ADP-ribose] polymerase 16) (PARP-16)
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type IV membrane protein .
Tissue Specificity	Brain,Carcinoma,
Function	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,similarity:Contains 1 PARP catalytic domain.,
Background	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,similarity:Contains 1 PARP catalytic domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PAR16 Monoclonal Antibody