



# RETNB Monoclonal Antibody

|                                    |                                                                                                                                                                                                                                                                                                                             |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-06996                                                                                                                                                                                                                                                                                                                |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                         |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                                                                                                                                             |
| <b>Applications</b>                | WB                                                                                                                                                                                                                                                                                                                          |
| <b>Gene Name</b>                   | RETNLB CCRG FIZZ2 HXCP2 RETNL2 UNQ408/PRO770                                                                                                                                                                                                                                                                                |
| <b>Immunogen</b>                   | Synthesized peptide derived from part region of human protein AA range: 1-50                                                                                                                                                                                                                                                |
| <b>Specificity</b>                 | RETNB Monoclonal Antibody detects endogenous levels of protein.                                                                                                                                                                                                                                                             |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                                                                                                                                                                                                              |
| <b>Source</b>                      | Monoclonal, Mouse,IgG                                                                                                                                                                                                                                                                                                       |
| <b>Purification</b>                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                        |
| <b>Dilution</b>                    | WB 1:500-1:2000                                                                                                                                                                                                                                                                                                             |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                        |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                |
| <b>Synonyms</b>                    | RETNLB CCRG;FIZZ1;FIZZ2;HXCP2;RELM beta;RELMb;RELMbeta;resisti; RETNB;RETNLB CCRG FIZZ2 HXCP2 RETNL2 UNQ408/PRO770;Resistin-like beta (Colon and small intestine-specific cysteine-rich protein) (Colon carcinoma-related gene protein) (Cysteine-rich secreted protein A12-alpha-like 1) (Cysteine-rich secreted protein F |
| <b>Observed Band</b>               | 12kD                                                                                                                                                                                                                                                                                                                        |
| <b>Calculated Molecular Weight</b> | 12kD                                                                                                                                                                                                                                                                                                                        |
| <b>Cell Pathway</b>                | Secreted.                                                                                                                                                                                                                                                                                                                   |
| <b>Tissue Specificity</b>          | Expressed only in the gastrointestinal tract, particularly the colon.                                                                                                                                                                                                                                                       |
| <b>Function</b>                    | function:Probable hormone.,similarity:Belongs to the resistin/FIZZ family., subunit:Homodimer; disulfide-linked.,tissue specificity:Expressed only in the gastrointestinal tract, particularly the colon.,                                                                                                                  |
| <b>Background</b>                  | function:Probable hormone.,similarity:Belongs to the resistin/FIZZ family., subunit:Homodimer; disulfide-linked.,tissue specificity:Expressed only in the gastrointestinal tract, particularly the colon.,                                                                                                                  |




**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**



Western Blot analysis of various cells using RETNB Monoclonal Antibody