



UR2R Monoclonal Antibody

Catalog No	YP-mAb-07558
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UTS2R GPR14
Immunogen	Synthesized peptide derived from human protein . at AA range: 130-210
Specificity	UR2R Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPR14;UR-2-R;UTR;UTR2;Senr;UR2R_HUMAN;UTS2R;G-protein coupled receptor 14;Urotensin II receptor (UR-II-R);UR2R_MOUSE;UR2R_RAT;G-protein coupled sensory epithelial neuropeptide-like receptor (SENr);urotensin 2 receptor;G protein-coupled receptor 14;UR2R;UTS2R GPR14;Urotensin-2 receptor (UR-2-R) (G-protein coupled receptor 14) (Urotensin II receptor) (UR-II-R)
Observed Band	42kD
Calculated Molecular Weight	42kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Most abundant expression in the heart and pancreas.
Function	function:High affinity receptor for urotensin-2 and urotensin-2B. The activity of this receptor is mediated by a G-protein that activate a phosphatidylinositol-calcium second messenger system.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Most abundant expression in the heart and pancreas.,



Background

function:High affinity receptor for urotensin-2 and urotensin-2B. The activity of this receptor is mediated by a G-protein that activate a phosphatidylinositol-calcium second messenger system.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Most abundant expression in the heart and pancreas.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using UR2R Monoclonal Antibody