




ZN273 mouse mAb

Catalog No	YP-mAb-08060
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF273
Immunogen	Synthesized peptide derived from human ZN273 AA range: 98-148
Specificity	This antibody detects endogenous levels of ZN273 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.175% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN273;ZNF273;Zinc finger protein 273 (Zinc finger protein HZF9)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Bladder,Bone marrow,Uterus,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 13 C2H2-type zinc fingers.,
Background	This gene is a member of the krueppel C2H2-type zinc-finger protein family and encodes a protein with 13 C2H2-type zinc fingers and a KRAB domain. This nuclear protein is involved in transcriptional regulation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZN273 mouse mAb