



ZN711 mouse mAb


Catalog No	YP-mAb-08179
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF711 CMPX1 ZNF6
Immunogen	Synthesized peptide derived from human ZN711 AA range: 438-488
Specificity	This antibody detects endogenous levels of ZN711 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.294% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CMPX1;dJ75N13.1;MRX97;Zfp711;Zinc finger protein 6 (CMPX1);Zinc finger protein 6;Zinc finger protein 711;ZN711;ZN711_HUMAN;ZNF4;ZNF5;ZNF6;ZNF711;ZNF711 CMPX1 ZNF6;Zinc finger protein 711 (Zinc finger protein 6)
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in neural tissues.
Function	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 11 C2H2-type zinc fingers.,
Background	This gene encodes a zinc finger protein of unknown function. It bears similarity to a zinc finger protein which acts as a transcriptional activator. This gene lies in a region of the X chromosome which has been associated with mental retardation. [provided by RefSeq, Jul 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZN711 mouse mAb