



## CYTSA mouse mAb

|                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-08505                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Applications</b>                | WB                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Gene Name</b>                   | SPECC1L CYTSA KIAA0376                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Immunogen</b>                   | Synthesized peptide derived from human CYTSA AA range: 150-200                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Specificity</b>                 | This antibody detects endogenous levels of CYTSA at Human/Mouse/Rat                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Source</b>                      | Monoclonal, Mouse,IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Purification</b>                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                                                                                                                          |
| <b>Dilution</b>                    | WB 1:500-1:2000                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Synonyms</b>                    | cytopsin A;CYTSA;KIAA0376;SPECC1 like protein;SPECC1L;SPECC1L CYTSA KIAA0376                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Observed Band</b>               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Calculated Molecular Weight</b> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Cell Pathway</b>                | Cytoplasm, cytoskeleton . Cytoplasm, cytoskeleton, spindle . Cell junction, gap junction . Colocalizes with acetylated alpha-tubulin, gamma-tubulin and F-actin. Also observed in a ring around gamma-tubulin containing centrioles possibly in the microtubule organizing center.                                                                                                                                                                                            |
| <b>Tissue Specificity</b>          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Function</b>                    | function:Involved in cytokinesis and spindle organization.,similarity:Belongs to the cytopsin-A family.,similarity:Contains 1 CH (calponin-homology) domain.,                                                                                                                                                                                                                                                                                                                 |
| <b>Background</b>                  | This gene encodes a coiled-coil domain containing protein. The encoded protein may play a critical role in actin-cytoskeletal reorganization during facial morphogenesis. Mutations in this gene are a cause of oblique facial clefting-1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. A read-through transcript composed of SPECC1L (sperm antigen with calponin homology and coiled-coil domains 1-like) and the |



downstream ADORA2A (adenosine A2a receptor) gene sequence has been identified, but it is thought to be non-coding. [provided by RefSeq, Jun 2013],


**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**



Western Blot analysis of various cells using CYTSA mouse mAb