



COQ2 mouse mAb

Catalog No	YP-mAb-08748
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	COQ2 CL640
Immunogen	Synthesized peptide derived from human COQ2 AA range: 160-210
Specificity	This antibody detects endogenous levels of COQ2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CL640;COQ2;COQ2 homolog;FLJ26072;hCOQ2;PHB:polyprenyltransferase; COQ2 CL640
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion inner membrane ; Multi-pass membrane protein ; Matrix side .
Tissue Specificity	Widely expressed. Present in all of the tissues tested. Expressed at higher level in skeletal muscle, adrenal glands and the heart.
Function	disease:Defects in COQ2 are a cause of coenzyme Q10 deficiency [MIM:607426]. Coenzyme Q10 deficiency is an autosomal recessive disorder with variable manifestations. It can be associated with three main clinical phenotypes: a predominantly myopathic form with central nervous system involvement, an infantile encephalomyopathy with renal dysfunction and an ataxic form with cerebellar atrophy.,function:Catalyzes the prenylation of para-hydroxybenzoate (PHB) with an all-trans polyprenyl group. Mediates the second step in the final reaction sequence of coenzyme Q (CoQ) biosynthesis, which is the condensation of the polyisoprenoid side chain with PHB.,pathway:Cofactor biosynthesis; ubiquinone biosynthesis.,similarity:Belongs to the ubiA prenyltransferase family.,tissue specificity:Widely expressed. Present in all of



the tissues tested. Expressed at higher level in skeletal muscle, adrenal gland

Background

This gene encodes an enzyme that functions in the final steps in the biosynthesis of CoQ (ubiquinone), a redox carrier in the mitochondrial respiratory chain and a lipid-soluble antioxidant. This enzyme, which is part of the coenzyme Q10 pathway, catalyzes the prenylation of parahydroxybenzoate with an all-trans polyprenyl group. Mutations in this gene cause coenzyme Q10 deficiency, a mitochondrial encephalomyopathy, and also COQ2 nephropathy, an inherited form of mitochondriopathy with primary renal involvement. [provided by RefSeq, Oct 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using COQ2 mouse mAb