



CHCH3 mouse mAb

Catalog No	YP-mAb-09824
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CHCHD3 MINOS3
Immunogen	Synthesized peptide derived from human CHCH3 AA range: 144-194
Specificity	This antibody detects endogenous levels of human CHCH3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CHCH3_HUMAN;CHCHD 3;Chchd3;Coiled coil helix coiled coil helix domain containing 3;Coiled coil helix coiled coil helix domain containing protein 3 mitochondrial;Coiled-coil-helix-coiled-coil-helix domain-containing protein 3; FLJ242;MINOS3;mitochondrial;Mitochondrial inner membrane organizing system 3;PPP1R22;Protein phosphatase 1 regulatory subunit 22;CHCH3; CHCHD3 MINOS3; mitochondrial
Observed Band	
Calculated Molecular Weight	26 kDa, 2 kDa; 26kD
Cell Pathway	Mitochondrion inner membrane ; Lipid-anchor ; Intermembrane side . Cytoplasm . Nucleus . Mitochondrion .
Tissue Specificity	Detected at low levels in brain, placenta, lung, liver, kidney and pancreas with increased levels in heart and skeletal muscle. Higher expression in primary lung cancers than in normal lung tissue.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images