



SHOX2 mouse mAb

Catalog No	YP-mAb-10999
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SHOX2 OG12X SHOT
Immunogen	Synthesized peptide derived from human SHOX2 AA range: 37-87
Specificity	This antibody detects endogenous levels of SHOX2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SHOT;Paired-related homeobox protein SHOT;OGI2X;OG12X;OG12;SHOX2;SHOX2 OG12X SHOT
Observed Band	
Calculated Molecular Weight	36kD
Cell Pathway	Nucleus.
Tissue Specificity	Expressed in heart, skeletal muscle, liver, lung, bone marrow fibroblast, pancreas and placenta.
Function	developmental stage:Expressed during cranofacial development as well as in heart.,function:May be a growth regulator and have a role in specifying neural systems involved in processing somatosensory information, as well as in face and body structure formation.,similarity:Belongs to the paired homeobox family. Bicoid subfamily.,similarity:Contains 1 homeobox DNA-binding domain.,similarity:Contains 1 OAR domain.,tissue specificity:Expressed in heart, skeletal muscle, liver, lung, bone marrow fibroblast, pancreas and placenta.,
Background	This gene is a member of the homeobox family of genes that encode proteins containing a 60-amino acid residue motif that represents a DNA binding domain. Homeobox genes have been characterized extensively as transcriptional



regulators involved in pattern formation in both invertebrate and vertebrate species. Several human genetic disorders are caused by aberrations in human homeobox genes. This locus represents a pseudoautosomal homeobox gene that is thought to be responsible for idiopathic short stature, and it is implicated in the short stature phenotype of Turner syndrome patients. This gene is considered to be a candidate gene for Cornelia de Lange syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using SHOX2 mouse mAb