




MCP1 mouse mAb

Catalog No	YP-mAb-11161
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MCP1, SCYA2,CCL2
Immunogen	Synthesized peptide derived from human MCP1. AA 50-64
Specificity	This antibody detects endogenous levels of MCP1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ccl2;MCP1;MCP-1;C-C motif chemokine 2;CCL 2;SCYA2;Monocyte chemoattractant protein 1;Monocyte chemotactic protein 1;Monocyte secretory protein JE; SCYA2
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MCP1 mouse mAb