



# PEF1 mouse mAb


|                                    |                                                                                                                                                                                                                                                                                                                                                                                                        |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-11202                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Applications</b>                | WB                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Gene Name</b>                   | PEF1 ABP32 UNQ1845/PRO3573                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Immunogen</b>                   | Synthesized peptide derived from human PEF1 AA range: 109-159                                                                                                                                                                                                                                                                                                                                          |
| <b>Specificity</b>                 | This antibody detects endogenous levels of PEF1 at Human/Mouse/Rat                                                                                                                                                                                                                                                                                                                                     |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                                |
| <b>Source</b>                      | Monoclonal, Mouse,IgG                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Purification</b>                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                                                   |
| <b>Dilution</b>                    | WB 1:500-1:2000                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Synonyms</b>                    | ABP32;PEF1;PEF1A;PEFLIN;PEF1 ABP32 UNQ1845/PRO3573                                                                                                                                                                                                                                                                                                                                                     |
| <b>Observed Band</b>               |                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Calculated Molecular Weight</b> |                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Cell Pathway</b>                | Cytoplasm . Endoplasmic reticulum . Membrane ; Peripheral membrane protein . Cytoplasmic vesicle, COPII-coated vesicle membrane ; Peripheral membrane protein . Membrane-associated in the presence of Ca(2+) (PubMed:11278427). Localizes to endoplasmic reticulum exit site (ERES) (By similarity). .                                                                                                |
| <b>Tissue Specificity</b>          |                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Function</b>                    | similarity:Contains 5 EF-hand domains.,subcellular location:Membrane-associated in the presence of Ca(2+).,subunit:Heterodimer; heterodimerizes with PDCD6/ALG2. Dissociates from PDCD6/ALG2 in presence of Ca(2+).,                                                                                                                                                                                   |
| <b>Background</b>                  | This gene encodes a calcium-binding protein belonging to the penta-EF-hand protein family. The encoded protein has been shown to form a heterodimer with the programmed cell death 6 gene product and may modulate its function in Ca(2+) signaling. Alternative splicing results in multiple transcript variants and a pseudogene has been identified on chromosome 1.[provided by RefSeq, May 2010], |

**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**

Western Blot analysis of various cells using PEF1 mouse mAb