




FOXD2 mouse mAb

Catalog No	YP-mAb-11364
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXD2 FKHL17 FREAC9
Immunogen	Synthesized peptide derived from human FOXD2 AA range: 177-227
Specificity	This antibody detects endogenous levels of FOXD2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FKHL17;Forkhead box protein D2;Forkhead-related protein FKHL17;Forkhead-related transcription;FOXD2;FOXD2 FKHL17 FREAC9
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Kidney specific.
Function	function:Probable transcription factor involved in embryogenesis and somatogenesis.,similarity:Contains 1 fork-head DNA-binding domain.,tissue specificity:Kidney specific.,
Background	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using FOXD2 mouse mAb