




ERLN1 mouse mAb

Catalog No	YP-mAb-12104
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ERLIN1 C10orf69 KE04 KEO4 SPFH1
Immunogen	Synthesized peptide derived from human ERLN1 AA range: 185-235
Specificity	This antibody detects endogenous levels of ERLN1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KE04;Erlin-1;Erlin 1;ER lipid raft associated 1;Endoplasmic reticulum lipid raft-associated protein 1;ERLN1;ERLIN1 C10orf69 KE04 KEO4 SPFH1
Observed Band	
Calculated Molecular Weight	38kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type II membrane protein . Associated with lipid raft-like domains of the endoplasmic reticulum membrane.
Tissue Specificity	Expressed in heart, placenta, liver, kidney, pancreas, prostate, testis, ovary and small intestine.
Function	similarity:Belongs to the band 7/mec-2 family.,subcellular location:Associated with lipid raft-like domains of the endoplasmic reticulum membrane.,tissue specificity:Expressed in heart, placenta, liver, kidney, pancreas, prostate, testis, ovary and small intestine.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ERLN1 mouse mAb