




PRR12 mouse mAb

Catalog No	YP-mAb-12359
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRR12 KIAA1205
Immunogen	Synthesized peptide derived from human PRR12 AA range: 1079-1129
Specificity	This antibody detects endogenous levels of PRR12 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NOC;KIAA1205;PRR12;PRR12 KIAA1205
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Cell junction, synapse, postsynaptic density . Cell junction, synapse, synaptosome .
Tissue Specificity	
Function	similarity:Contains 2 A.T hook DNA-binding domains.,
Background	This gene encodes a proline-rich protein that contains two A-T hook DNA binding domains. A chromosomal translocation and gene fusion between this gene and zinc finger, MIZ-type containing 1 (Gene ID: 57178) may underlie intellectual disability and neuropsychiatric problems in a human patient. Enriched expression of this gene in embryonic mouse brain suggests that this gene may play a role in nervous system development. [provided by RefSeq, Jul 2016],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using PRR12 mouse mAb