




NCKP5 mouse mAb

Catalog No	YP-mAb-12405
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NCKAP5 ERIH NAP5
Immunogen	Synthesized peptide derived from human NCKP5 AA range: 965-1015
Specificity	This antibody detects endogenous levels of NCKP5 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NAP5;Peripheral clock protein;Nck associated protein 5;NAP 5;ERIH2;NCKP5; NCKAP5 ERIH NAP5
Observed Band	210kD
Calculated Molecular Weight	210kD
Cell Pathway	
Tissue Specificity	Expressed in fetal and adult brain, leukocytes and fetal fibroblasts.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NCKP5 mouse mAb