



Cables2 Monoclonal Antibody


Catalog No	YP-mAb-12692
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CABLES2
Immunogen	The antiserum was produced against synthesized peptide derived from human Ik3-2. AA range:91-140
Specificity	Cables2 Monoclonal Antibody detects endogenous levels of Cables2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2orf15;CABLES 2;Cdk5 and Abl enzyme substrate 2;CDK5 and ABL1 enzyme substrate 2;dJ98M14.2;Ik3 2;Interactor with CDK3 2;Cables2;C20orf150;Ik3-2
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	
Tissue Specificity	Eye,
Function	function:Unknown. Probably involved in G1-S cell cycle transition., similarity:Belongs to the cyclin family.,subunit:Binds to CDK3, CDK5 and ABL1. The C-terminal cyclin-box-like region binds to CDK5.,
Background	function:Unknown. Probably involved in G1-S cell cycle transition., similarity:Belongs to the cyclin family.,subunit:Binds to CDK3, CDK5 and ABL1. The C-terminal cyclin-box-like region binds to CDK5.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Cables2 Monoclonal Antibody