



NTR2 Monoclonal Antibody

Catalog No	YP-mAb-12782
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NTSR2
Immunogen	The antiserum was produced against synthesized peptide derived from human NTR2. AA range:151-200
Specificity	NTR2 Monoclonal Antibody detects endogenous levels of NTR2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NTR2;NTSR2;Neurotensin receptor type 2;NT-R-2;Levocabastine-sensitive neurotensin receptor
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Expressed in prostate (at protein level).
Function	function:Receptor for the tridecapeptide neurotensin. It is associated with G proteins that activate a phosphatidylinositol-calcium second messenger system., online information:Neurotensin receptor entry,similarity:Belongs to the G-protein coupled receptor 1 family.,
Background	The protein encoded by this gene belongs to the G protein-coupled receptor family that activate a phosphatidylinositol-calcium second messenger system. Binding and pharmacological studies demonstrate that this receptor binds neurotensin as well as several other ligands already described for neurotensin NT1 receptor. However, unlike NT1 receptor, this gene recognizes, with high affinity, levocabastine, a histamine H1 receptor antagonist previously shown to



compete with neurotensin for low-affinity binding sites in brain. These activities suggest that this receptor may be of physiological importance and that a natural agonist for the receptor may exist. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

