



SR-5A Monoclonal Antibody

Catalog No	YP-mAb-12820
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HTR5A
Immunogen	The antiserum was produced against synthesized peptide derived from human HTR5A. AA range:15-64
Specificity	SR-5A Monoclonal Antibody detects endogenous levels of SR-5A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	5 HT 5;5 HT 5A;5 HT5A;HTR5A;Serotonin receptor 5A;SR-5A;5-hydroxytryptamine receptor 5A;5-HT-5;5-HT-5A;5-HT5A
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Brain,Cerebellum,
Function	function:This is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. The activity of this receptor is mediated by G proteins., similarity:Belongs to the G-protein coupled receptor 1 family.,
Background	The neurotransmitter serotonin (5-hydroxytryptamine, 5-HT) has been implicated in a wide range of psychiatric conditions and also has vasoconstrictive and vasodilatory effects. The gene described in this record is a member of 5-hydroxytryptamine (serotonin) receptor family and encodes a multi-pass membrane protein that functions as a receptor for 5-hydroxytryptamine and couples to G-proteins. This protein has been shown to function in part through



the regulation of intracellular Ca²⁺ mobilization. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SR-5A Monoclonal Antibody

