



FPR1 Monoclonal Antibody

Catalog No	YP-mAb-13250
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FPR1
Immunogen	The antiserum was produced against synthesized peptide derived from human FPR1. AA range:151-200
Specificity	FPR1 Monoclonal Antibody detects endogenous levels of FPR1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	fMet Leu Phe receptor;FMLP;fMLP receptor;formyl peptide receptor 1;FPR; FPR1;fMet-Leu-Phe receptor;N-formyl peptide receptor;N-formylpeptide chemoattractant receptor
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Internalizes in presence of its ligands, fMLP, TAFA4 and CTSG. .
Tissue Specificity	Neutrophils.
Function	function:High affinity receptor for N-formyl-methionyl peptides, which are powerful neutrophils chemotactic factors. Binding of FMLP to the receptor causes activation of neutrophils. This response is mediated via a G-protein that activates a phosphatidylinositol-calcium second messenger system., PTM:Phosphorylated; which is necessary for desensitization.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Neutrophils.,



Background

formyl peptide receptor 1(FPR1) Homo sapiens This gene encodes a G protein-coupled receptor of mammalian phagocytic cells that is a member of the G-protein coupled receptor 1 family. The protein mediates the response of phagocytic cells to invasion of the host by microorganisms and is important in host defense and inflammation.[provided by RefSeq, Jul 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using FPR1 Monoclonal Antibody