



GPR160 Monoclonal Antibody

Catalog No	YP-mAb-13314
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPR160
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR160. AA range:289-338
Specificity	GPR160 Monoclonal Antibody detects endogenous levels of GPR160 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPR160;Probable G-protein coupled receptor 160;GPCR150;G-protein coupled receptor GPCR1;hGPCR1
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Cervix,Pancreas,
Function	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

