



GPR31 Monoclonal Antibody

Catalog No	YP-mAb-13330
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPR31
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR31. AA range:181-230
Specificity	GPR31 Monoclonal Antibody detects endogenous levels of GPR31 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPR31;12-(S)-hydroxy-5;8;10;14-eicosatetraenoic acid receptor;12-(S)-HETE receptor;12-HETER;G-protein coupled receptor 31
Observed Band	33kD
Calculated Molecular Weight	33kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	G-protein coupled receptors,Placenta,
Function	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using GPR31
Monoclonal Antibody

