



GPR32 Monoclonal Antibody

Catalog No	YP-mAb-13331
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPR32
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR32. AA range:151-200
Specificity	GPR32 Monoclonal Antibody detects endogenous levels of GPR32 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DRV1;RVDR1;GPR32;Probable G-protein coupled receptor 32
Observed Band	33kD
Calculated Molecular Weight	33kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	
Function	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,
Background	This gene is intronless and encodes a member of the G-protein coupled receptor 1 family. The encoded protein binds to resolvins D1 and lipoxin A4 and has been linked to pulmonary inflammation. A related pseudogene has been identified on chromosome 19. [provided by RefSeq, Nov 2012],
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

