




Mic-1 Monoclonal Antibody

Catalog No	YP-mAb-13945
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C18orf8
Immunogen	The antiserum was produced against synthesized peptide derived from human MIC1. AA range:591-640
Specificity	Mic-1 Monoclonal Antibody detects endogenous levels of Mic-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C18orf8;MIC1;Colon cancer-associated protein Mic1;HsT2591;Mic 1;Mic-1; Uncharacterized protein C18orf8
Observed Band	75kD
Calculated Molecular Weight	75kD
Cell Pathway	Lysosome membrane . Late endosome membrane .
Tissue Specificity	Colon carcinoma,Epithelium,Lung,Synovium,
Function	similarity:Contains 1 Mic1 domain.,
Background	This gene encodes a colon cancer associated protein. [provided by RefSeq, Jan 2013],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Mic-1
Monoclonal Antibody