



CD66c/d Monoclonal Antibody

Catalog No	YP-mAb-14057
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CEACAM3/CEACAM6
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human CEACAM3/CEACAM6. AA range:31-80
Specificity	CD66c/d Monoclonal Antibody detects endogenous levels of CD66c/d protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD66c/d;CEACAM3/CEACAM6;Carcinoembryonic antigen-related cell adhesion molecule 3/Carcinoembryonic antigen-related cell adhesion molecule 6; CEACAM3;CD66D;CGM1;Carcinoembryonic antigen-related cell adhesion molecule 3;Carcinoembryonic antigen CGM1;CEACAM6;NCA; Carcinoembryonic antigen-related cell adhesion molecule 6;Non-specific crossreacting antigen;Normal cross-reacting antigen;CD66c
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Membrane; Single-pass type I membrane protein.
Tissue Specificity	CGM1a, the predominant CGM1 transcript, is granulocyte-specific. Not detected out of the granulocytic lineage, such as monocytes, lymphocytes, spleen, testis, colon, brain, liver, pancreas, thymus, ovary, placenta, skeletal muscle, prostate, small intestine, heart, lung and kidney.
Function	caution:This is not the ortholog of rat CEACAM3.,domain:The cytosolic domain is involved in S100A9 interaction.,function:Major granulocyte receptor mediating recognition and efficient opsonin-independent phagocytosis of CEACAM-binding microorganisms, including Neissiria, Moxarella and Haemophilus species, thus



playing an important role in the clearance of pathogens by the innate immune system. Responsible for RAC1 stimulation in the course of pathogen phagocytosis.,PTM:Tyrosine-phosphorylated in response to microbial binding. Tyr-230 and Tyr-241 are both required for phosphorylation to be detected., similarity:Belongs to the immunoglobulin superfamily. CEA family., similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain., subunit:Interacts with S100A9/calprotectin. This interaction is calcium-dependent, but independent of CEACAM3 phosphorylation.,tissue specificity:CGM1a, the pred

Background

This gene encodes a member of the family of carcinoembryonic antigen-related cell adhesion molecules (CEACAMs), which are used by several bacterial pathogens to bind and invade host cells. The encoded transmembrane protein directs phagocytosis of several bacterial species that is dependent on the small GTPase Rac. It is thought to serve an important role in controlling human-specific pathogens by the innate immune system. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2013],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using CD66c/d Monoclonal Antibody

