



HXK I Monoclonal Antibody

Catalog No	YP-mAb-14778
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HK1
Immunogen	The antiserum was produced against synthesized peptide derived from human HXK1. AA range:31-80
Specificity	HXK I Monoclonal Antibody detects endogenous levels of HXK I protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BB4413;Brain form hexokinase;dea;DrHXK1;EC 2.7.1.1;Glycolytic enzyme; HEXOKIN;hexokinase I;Hexokinase PI;Hexokinase type I;Hexokinase;tumor isozyme;Hexokinase-1;Hexokinase-A;HK I;HK1;HK1 tb;Hk1-s;HK1-ta;HK1-tb; HK1-tc;HKD;HKI;HMSNR;HXK1;HXK1_HUMAN;im:7148527;mHk1-s;wu:fc9d8; wu:fc16e2;wu:fc21e2;wu:fc14b11;zgc:5579;zgc:77618;HXK I
Observed Band	109kD
Calculated Molecular Weight	109kD
Cell Pathway	Mitochondrion outer membrane ; Peripheral membrane protein . Cytoplasm, cytosol . The mitochondrial-binding peptide (MBP) region promotes association with the mitochondrial outer membrane (Probable). Dissociates from the mitochondrial outer membrane following inhibition by N-acetyl-D-glucosamine, leading to relocation to the cytosol (PubMed:27374331). .
Tissue Specificity	Isoform 2: Erythrocyte specific (Ref.6). Isoform 3: Testis-specific (PubMed:10978502). Isoform 4: Testis-specific (PubMed:10978502).
Function	catalytic activity:ATP + D-hexose = ADP + D-hexose 6-phosphate., disease:Defects in HK1 are the cause of hexokinase deficiency [MIM:235700]. Hexokinase deficiency is a rare autosomal recessive disease with nonspherocytic hemolytic anemia as the predominant clinical feature.,



domain: The N- and C-terminal halves of this hexokinase show extensive sequence similarity to each other. The catalytic activity is associated with the C-terminus while regulatory function is associated with the N-terminus., enzyme regulation: Hexokinase is an allosteric enzyme inhibited by its product Glc-6-P., miscellaneous: In vertebrates there are four major glucose-phosphorylating isoenzymes, designated hexokinase I, II, III and IV (glucokinase)., online information: Hexokinase entry, pathway: Carbohydrate metabolism; hexose metabolism., similarity: Belongs to the hexokinase family., subcellular location: Its hydrophobic N-ter

Background

Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes a ubiquitous form of hexokinase which localizes to the outer membrane of mitochondria. Mutations in this gene have been associated with hemolytic anemia due to hexokinase deficiency. Alternative splicing of this gene results in several transcript variants which encode different isoforms, some of which are tissue-specific. [provided by RefSeq, Apr 2016],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using HXK I Monoclonal Antibody