



NIM1 Monoclonal Antibody

Catalog No	YP-mAb-14879
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NIM1
Immunogen	The antiserum was produced against synthesized peptide derived from human NIM1. AA range:347-396
Specificity	NIM1 Monoclonal Antibody detects endogenous levels of NIM1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Hypothetical protein LOC167359;Hypothetical protein MGC4215;Nim1; NIM1_HUMAN;Serine/threonine-protein kinase NIM1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,cytoplasm,integral component of membrane,
Tissue Specificity	Brain,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein., cofactor:Magnesium.,enzyme regulation:Activated by phosphorylation at Thr-229, probably by autophosphorylation.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein., cofactor:Magnesium.,enzyme regulation:Activated by phosphorylation at Thr-229, probably by autophosphorylation.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using NIM1 Monoclonal Antibody