



PP2C κ Monoclonal Antibody

Catalog No	YP-mAb-14953
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PPM1K
Immunogen	The antiserum was produced against synthesized peptide derived from human PPM1K. AA range:205-254
Specificity	PP2C κ Monoclonal Antibody detects endogenous levels of PP2C κ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PP2C κ ;PP2C μ ;PPM1K;PTMP;UG0882E07;PP2C κ ;Protein phosphatase 1K mitochondrial;Protein phosphatase 1K;mitochondrial;PP2C domain-containing protein phosphatase 1K;PP2C-like mitochondrial protein;PP2C-type mitochondrial phosphoprotein phosphatase;Protein phosphatase 2C isoform κ ;PP2C- κ
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	Amygdala,Cerebellum,Lymph node,Pancreas,Testis,
Function	catalytic activity:A phosphoprotein + H(2)O = a protein + phosphate., caution:Ref.9 has crystallized PPM1K in the presence of magnesium ions. However, PubMed:17336929 reported that no activity toward p-nitrophenylphosphate was seen in the absence of manganese ions and magnesium could not substitute for manganese.,cofactor: Binds 1 magnesium or manganese ion per subunit.,function:Regulates the mitochondrial permeability transition pore and is essential for cellular survival and development., similarity:Belongs to the PP2C family.,similarity:Contains 1 PP2C-like domain.,



Background

This gene encodes a member of the PPM family of Mn^{2+}/Mg^{2+} -dependent protein phosphatases. The encoded protein, essential for cell survival and development, is targeted to the mitochondria where it plays a key role in regulation of the mitochondrial permeability transition pore. [provided by RefSeq, Sep 2012],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using PP2C α Monoclonal Antibody