



PRKY Monoclonal Antibody

Catalog No	YP-mAb-14961
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PRKY
Immunogen	The antiserum was produced against synthesized peptide derived from human PRKY. AA range:61-110
Specificity	PRKY Monoclonal Antibody detects endogenous levels of PRKY protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRKY;protein kinase;Y linked;Serine/threonine-protein kinase PRKY
Observed Band	34kD
Calculated Molecular Weight	34kD
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	
Background	This gene is similar to the protein kinase, X-linked gene in the pseudoautosomal region of the X chromosome. The gene is classified as a transcribed pseudogene because it has lost a coding exon that results in all transcripts being candidates for nonsense-mediated decay (NMD) and unlikely to express a protein. Abnormal recombination between this gene and a related gene on chromosome X is a frequent cause of XX males and XY females. [provided by RefSeq, Jul 2010],
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

