



GHITM Monoclonal Antibody


Catalog No	YP-mAb-15917
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GHITM
Immunogen	The antiserum was produced against synthesized peptide derived from human GHITM. AA range:31-80
Specificity	GHITM Monoclonal Antibody detects endogenous levels of GHITM protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Dermal papilla derived protein 2;Dermal papilla-derived protein 2;DERP2;Ghitm; GHITM_HUMAN;Growth hormone inducible transmembrane protein;Growth hormone-inducible transmembrane protein;HSPC282;MICS1;Mitochondrial morphology and cristae structure 1;My21;PTD1;TMBIM5;Transmembrane BAX inhibitor motif containing 5;Transmembrane BAX inhibitor motif containing protein 5;Transmembrane BAX inhibitor motif-containing protein 5;My021
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Mitochondrion inner membrane ; Multi-pass membrane protein .
Tissue Specificity	Blood,Brain,Fetal brain,Hair follicle dermal papilla,Kidney
Function	similarity:Belongs to the BI1 family.,
Background	similarity:Belongs to the BI1 family.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using GHITM Monoclonal Antibody