



KV3.2 Monoclonal Antibody

Catalog No	YP-mAb-16454
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KCNC2
Immunogen	The antiserum was produced against synthesized peptide derived from human Potassium Channel Kv3.2b. AA range:589-638
Specificity	KV3.2 Monoclonal Antibody detects endogenous levels of KV3.2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KCNC2;KSHIIIA;Kv 3.2;Kv3.2d voltage gated potassium channel;Potassium channel Kv3.2;Shaw related voltage gated potassium channel protein 2;KV3.2; Potassium voltage-gated channel subfamily C member 2;Voltage-gated potassium channel Kv3.2
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Membrane ; Multi-pass membrane protein . Perikaryon . Cell projection, axon . Cell projection, dendrite . Cell junction, synapse, postsynaptic cell membrane . Cell junction, synapse, presynaptic cell membrane . Cell junction, synapse, synaptosome . Cell junction, synapse . Apical cell membrane . Basolateral cell membrane . Colocalizes with parvalbumin in globus pallidus neurons. Localizes in thalamocortical axons and synapses. Localizes on the surface of cell somata, proximal dendrites and axonal membranes. Also detected throughout the neuropil. Localized in starburst cell somata and proximal dendrite processes. Colocalized with GABA in presynaptic terminals. Clustered in patches in somatic and proximal dendritic membrane as well as in axons an
Tissue Specificity	Amygdala,Brain,



Function

domain:The segment S4 is probably the voltage-sensor and is characterized by a series of positively charged amino acids at every third position.,domain:The tail may be important in modulation of channel activity and/or targeting of the channel to specific subcellular compartments.,function:Mediates the voltage-dependent potassium ion permeability of excitable membranes. Assuming opened or closed conformations in response to the voltage difference across the membrane, the protein forms a potassium-selective channel through which potassium ions may pass in accordance with their electrochemical gradient. Channel properties are modulated by subunit assembly.,similarity:Belongs to the potassium channel family. C (Shaw) subfamily.,subunit:The potassium channel is probably composed of a homo- or heterotetrameric complex of pore-forming subunits that can associate with modulating accessory subun

Background

The Shaker gene family of Drosophila encodes components of voltage-gated potassium channels and is comprised of four subfamilies. Based on sequence similarity, this gene is similar to one of these subfamilies, namely the Shaw subfamily. The protein encoded by this gene belongs to the delayed rectifier class of channel proteins and is an integral membrane protein that mediates the voltage-dependent potassium ion permeability of excitable membranes. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using KV3.2 Monoclonal Antibody