



PCTAIRE-2 Monoclonal Antibody

Catalog No	YP-mAb-16767
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CDK17
Immunogen	Synthesized peptide derived from PCTAIRE-2 . at AA range: 70-150
Specificity	PCTAIRE-2 Monoclonal Antibody detects endogenous levels of PCTAIRE-2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	643598J1Rik;CDK17;CDK17_HUMAN;Cell division protein kinase 17;Cyclin-dependent kinase 17;EC 2.7.1.37;MGC2519;PCTAIRE-motif protein kinase 2;PCTK2;Protein kinase cdc2-related PCTAIRE-2;Serine/threonine-protein kinase PCTAIRE-2;PCTAIRE-2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Hippocampus,
Function	similarity:Belongs to the protein kinase superfamily.,
Background	The protein encoded by this gene belongs to the cdc2/cdkx subfamily of the ser/thr family of protein kinases. It has similarity to a rat protein that is thought to play a role in terminally differentiated neurons. Alternatively spliced transcript variants encoding different isoforms have been found. [provided by RefSeq, Jul 2010],



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

