



Connexin 47 Monoclonal Antibody

Catalog No	YP-mAb-17006
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GJC2
Immunogen	The antiserum was produced against synthesized peptide derived from human CXG2. AA range:21-70
Specificity	Connexin 47 Monoclonal Antibody detects endogenous levels of Connexin 47 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	B23382L12Rik;Connexin 46.6;Connexin 47;Connexin-46.6;Connexin-47;Connexin46.6;Connexin47;CX 46.6;Cx 47;Cx46.6;Cx47;CXG2_HUMAN;Gap junction alpha 12 protein;Gap junction alpha-12 protein;Gap junction gamma 2 protein;Gap junction gamma-2 protein;Gap junction protein alpha 12 47kDa;Gap junction protein gamma 2 47kDa;GAP JUNCTION PROTEIN;47-KD;gamma 12;47kDa;gamma 2;GJA 12;GJA12;GJC 2;Gjc2;HLD 2;HLD2;LMPH1C;MGC15119;MGC1356;OTTHUMP37763;PMLDAR;SPG44
Observed Band	47kD
Calculated Molecular Weight	47kD
Cell Pathway	Cell membrane; Multi-pass membrane protein. Cell junction, gap junction.
Tissue Specificity	Expressed in central nervous system, in sciatic nerve and sural nerve. Also detected in skeletal muscles.
Function	caution:It is uncertain whether Met-1 or Met-4 is the initiator.,disease:Defects in GJC2 are the cause of Leukodystrophy hypomyelinating type 2 (HLD2) [MIM:608804]; also known as Pelizaeus-Merzbacher-like disease autosomal recessive type 1. HLD2 is an autosomal recessive hypomyelinating leukodystrophy characterized by nystagmus, impaired motor development,



ataxia, choreoathetotic movements, dysarthria and progressive spasticity., function:One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell. May play a role in myelination in central and peripheral nervous systems.,similarity:Belongs to the connexin family. Gamma-type subfamily.,subunit:A connexon is composed of a hexamer of connexins. Interacts with TJP1.,tissue specificity:Expressed in central nervous system,

Background

This gene encodes a gap junction protein. Gap junction proteins are members of a large family of homologous connexins and comprise 4 transmembrane, 2 extracellular, and 3 cytoplasmic domains. This gene plays a key role in central myelination and is involved in peripheral myelination in humans. Defects in this gene are the cause of autosomal recessive Pelizaeus-Merzbacher-like disease-1. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Connexin 47 Monoclonal Antibody